



# Thermo Scientific HeraSafe 2030i Biological Safety Cabinet

Setting new standards in biological safety cabinets





# Smarter protection

that never takes a day off

The Herasafe 2030i series of biological safety cabinets leverages a legacy of global leadership in BSC design with breakthrough innovations to optimize contamination prevention, user convenience and information management, delivering the most technologically advanced Class II cabinets available.

Equipped with “smart” self monitoring safety features, intuitive graphic touchscreen interface and built in connectivity to preserve important performance data, the HeraSafe 2030i delivers performance with unmatched precision, for demanding cell culture applications and other sensitive work requiring assured safety for product, personnel and environment, especially for bio pharma and advanced research initiatives.

Thermo Scientific biological safety cabinets are designed to maximize sample protection and user safety, with an emphasis on containment, comfort and convenience, all day, every day. For more than 5 decades we have been a pioneer and innovator of clean air technologies. Our global vision and experience have enabled us to continually provide the best application based solutions ideally suited for your requirements and working environment, with compliance to all major worldwide safety standards.



## Precise self adjusting protection

Autonomous safety features maintain proper air flow and monitor critical conditions in real time without user interaction. This supports valuable research enabling our customers to focus where needed.



## Simplified operation and ergonomics

Readily understandable graphical display provides the information you need, while electronic sash with fully hinged opening enables easy access for supplies and equipment, and simplifies easy cleaning and disinfection



## Enhanced information management

Electronically collect key operational and performance data on the secure Thermo Scientific Cloud and remotely monitor from anywhere, to evaluate your processes and accelerate your discoveries



# Simplified operation and ergonomics



Readily understandable graphical display provides the information you need, while electronic sash with fully hinged opening enables easy access for supplies and equipment, and enables easy cleaning and disinfection.

Helping ensure that the BSC is managing safety conditions autonomously, and is simple and intuitive to use allows our customers to fully focus on demanding applications and precision science.

## | Intuitive | Simple | Convenient

Real time display & BSC naming capability

Company common menu bar

Globally recognizable icons

Animated airflow

Company common "Health Status" icon

Additional velocity information

Frequent programs

Push button operations

Programmable functions:  
Edit auto start

- Acme Labs BSC 2
- 11:22:24AM | May 12, 2015
- Edit AutoStart
- Window Open
- Temporarily Mute Window Alarm
- Fan On
- Light On
- Receptacle On
- Save



### Smart Clean Plus

Fully opening front-hinged window design enables easy cleaning and disinfection with full access to the interior cabinet space.

Full opening also allows convenient posture for thorough cleaning and visible reach to all cabinet surfaces including rear wall.



### 316L Stainless steel interior (standard)

Most resistant to scratches and chemical attack in the user zone. All surfaces inside the HeraSafe 2030i are from 316L stainless steel. The finish is "brushed and sanded" which both reduces glare and offers better scratch resistance.

# Smart technology delivering optimal containment and energy efficiency

Autonomous safety features insure continuous accuracy that never takes a day off – providing optimal protection from contamination by compensating to match the way you work and lab environment

Sash Opening 200/250mm | Smart Flow+ and Dave+ | Password Protection  
| UV Protection Energy | Efficient Design |



## Electronic front sash window


with push button simplicity the front window automatically positions itself at the position mandated by local regulation to assure proper working conditions. Sloped front window allows fatigue free working position and is aerosol tight sealed when closed

**200mm sash opening (Standard)**  
provides optimum balance of comfort plus efficiency.

**250mm sash opening (optional)**  
provides additional convenience for quick and safe consumables movement.



## Class II BSC



## SmartFlow Plus

patented proprietary design automatically compensates unique energy saving dual DC blower motors to insure the proper airflow rates no matter the loading condition of the protective HEPA filters. Unlike typical BSC's which rarely perform at the same levels day after day due to HEPA loading, and potential disruptions to airflow such as paper covering vent openings.

### How does this work?

Independent intelligent DC motors for downflow and inflow monitor and control fan speed in real time to maintain balanced airflow during filter loading or temporary obstruction.

Exclusive to our BSCs, independent supply and exhaust blowers automatically balance both the downflow and inflow velocities to maintain both operator and sample protection.

## Dave Plus

our unique digital airflow verification system validates protection by monitoring airflow velocities with independent pressure sensors which detect and alert to any changes across the exhaust and downflow plenums.

### How does this work?

By sensing pressure in the throat of the blowers, DAVEPlus monitors airflow directly for the inflow and downflow in the Herasafe 2030i to assure both personal and sample protection.




## Controlled Passcode Access

for convenience the BSC can be configured via the user interface to operate in “open mode” or “secure mode” depending upon working preferences and lab environment. In secure mode a selectable passcode must be entered to enable operation of the BSC. In this way different units utilized by different personnel for specific applications, may be controlled by providing a unique passcode to each, or the same code to all.



## Self adjusting germicidal UV light intensity

for those who routinely employ UV germicidal light to their work surface at the end of the work day, we provide a time adjustable functionality to auto shut off the light at the desired interval to preserve bulb life. Our exclusive self adjusting exposure capability insures reproducible germicidal effectiveness with every use, by assuring performance at consistent levels to protect your important work.




## Energy efficient design

- Energy conserving dual DC motor design, generates less heat than older AC designs with double the usage lifetime and digitally controllable for precision performance
- Stand-by mode automatically reduces blower speed to protect samples with reduced energy usage until you are ready to resume work

### Optimizing total cost of operation

Thermo Scientific SmartFlow airflow compensation and balanced filters have led to safer BSCs that are easier to use and significantly less costly to operate.\*\*

### The more you use it, the more you save!



Costs, published in American Biotechnology Laboratory, September 2008.

\*\* Compared to previous models.


\*\*\* Comparison based on 15 year life of 1.2 m wide biological safety cabinets and 2011 average costs of energy and cooling in the EuroZone. "Traditional BSC" estimates based on published product specifications for readily available BSCs using AC motor technology.

# Innovative productivity enhancements simplify use, comfort and cleaning



## Full color touch screen user interface – safety you can see

Our advanced yet simple to operate, best in class graphic user interface is customer tested to provide easy navigation and convenient illustration of airflow safety and usage data to address any lingering user concerns about the cabinet functioning properly to protect them and their valuable samples. Performance details such as user interaction logs, errors and event logs are easily identified to drive continuous process improvement and reproducibility of results.



### Standard smart ports

provided as standard with every cabinet, Smart Ports provide clean and safe approach to incorporating vacuum tubing and cables through the side wall of the BSC, for improved organization and work efficiency. The side wall plenum is under negative pressure preventing possibility of leakage in or out of the BSC during operation.

### Customisation

options are available for the Herasafe 2030i. We have vast experience in designing custom adaptations to our BSC's. We can design from scratch or provide frequently requested solutions.

### Large interior work space/premium 316 non-corrosive stainless steel

conveniently large interior workspace provides ample room for incorporation of instruments and process optimization devices including robotic arms. We utilize a “brushed and sanded” finish to provide maximum scratch resistance and non-reflective properties.

### Padded armrests

adjustable ergonomic armrests are provided with each unit to assure proper arm support to reduce strain during prolonged work periods while extending reach within the cabinet. Removable gel filled comfort pads are available as an option to enhance the user experience.

### Interior LED lights

visibility within the cabinet is enhanced with the use of bright, energy efficient lighting, eliminating the potential for fatigue often experienced with use of standard fluorescent lighting

### Exclusive one touch start /stop function

exclusive one touch auto start and stop function activates or deactivate all cabinet components efficiently, and ready to use without lengthy wait times.

### Smart Clean hinged front window

our exclusive fully opening front hinged window design enables easy cleaning and full access to the cabinet interior work space for thorough cleaning/disinfection and incorporation of ancillary work tools and equipment as may be desired.





## Enhanced information management with the convenience of connectivity



Most of the enhanced information provided by the touch screen interface, can be remotely monitored and captured by web based or mobile app for subsequent analysis, eliminating the need for paper notebooks

**Standard USB data outputs** – enable export of data to external data systems if desired. A powered connection is included in the front panel to enable convenient accessory charging and accommodate future software update uploads

**Convenient event and alert messaging logs** – accessible in the user interface provide useful and important usage data and alert information for process optimization and maintenance

**Embedded connectivity** – ready for plug and play into the Thermo Fisher secure cloud as desired. Enables simplified information sharing with colleagues and enhances productivity through performance monitoring from anywhere, anytime improving process control, reporting and asset management, with use of the exclusive free Thermo Scientific mobile app

**Future proof** – once you have incorporated the Thermo Scientific Cloud you are positioned to benefit from continuous improvement of data management, as your capabilities with expand as the Thermo Scientific cloud evolves to include new functionality over time (active Q4 2018)



| Herasafe 2030i Ordering Information   |          |          | EN 12469                      |             |             |             | NSF/ANSI 49                 |             |
|---|----------|----------|-------------------------------|-------------|-------------|-------------|-----------------------------|-------------|
| Size/Internal width   | m (ft)   |          | 0.9 (3)                       | 1.2 (4)     | 1.5 (5)     | 1.8 (6)     | 1.2 (4)                     | 1.8 (6)     |
| 200mm Sash with Cross Beam UV-C in side walls   |          |          | 51032328                      | 51032329    | 51032330    | 51032331    |                             |             |
| 200mm Sash – no options included**<br>(250mm/10" sash can be optioned)  |          |          | 51032324                      | 51032325    | 51032326    | 51032327    |                             |             |
| 250mm (10") Sash with Cross Beam UV-C in side walls, indented single piece work tray and prep for 3x rear wall mount service taps |          |          |                               |             |             |             | 51032334                    | 51032335    |
| 250mm (10") Sash, indented single piece work tray   |          |          |                               |             |             |             | 51032332                    | 51032333    |
| Window operating aperture 10" (250mm)   | Option   |          | 51901247**                    |             |             |             | Standard                    |             |
| <b>Dimensions</b>   |          |          |                               |             |             |             |                             |             |
| External Dimension (Cabinet only)   | depth    | mm (in.) | 800 (31.5)                    | 800 (31.5)  | 800 (31.5)  | 800 (31.5)  | 800 (31.5)                  | 800 (31.5)  |
|   | width    | mm (in.) | 1000 (39.4)                   | 1300 (51.2) | 1600 (63.0) | 1900 (74.8) | 1300 (51.2)                 | 1900 (74.8) |
|   | height   | mm (in.) | 1536 (60.5)                   | 1536 (60.5) | 1536 (60.5) | 1536 (60.5) | 1536 (60.5)                 | 1536 (60.5) |
| Usable Work Surface<br>(net, without air intake grill)  | depth    | mm (in.) | 465 (18.3)                    | 465 (18.3)  | 465 (18.3)  | 465 (18.3)  | 465 (18.3)                  | 465 (18.3)  |
|   | width    | mm (in.) | 900 (35.4)                    | 1200 (47.2) | 1500 (59.1) | 1800 (70.9) | 1200 (47.2)                 | 1800 (70.9) |
| Shipping Dimension  | depth    | mm (in.) | 925 (36.4)                    | 925 (36.4)  | 925 (36.4)  | 925 (36.4)  | 925 (36.4)                  | 925 (36.4)  |
|   | width    | mm (in.) | 1110 (43.7)                   | 1410 (55.5) | 1710 (67.3) | 2010 (79.1) | 1410 (55.5)                 | 2010 (79.1) |
|   | height   | mm (in.) | 1710 (67.3)                   | 1710 (67.3) | 1710 (67.3) | 1710 (67.3) | 1710 (67.3)                 | 1710 (67.3) |
| Working height of front window  | mm (in.) |          | 200 (8)                       |             |             |             | 250 (10)                    |             |
| Maximum window opening<br>(for cleaning, loading/unloading)   | mm (in.) |          | 773 (30.4)                    | 773 (30.4)  | 773 (30.4)  | 773 (30.4)  | 773 (30.4)                  | 773 (30.4)  |
| Work surface height   | mm (in.) |          | 750-950 (29.5-37.4)           |             |             |             | 750-950 (29.5-37.4)         |             |
| <b>Weight</b>   |          |          |                               |             |             |             |                             |             |
| Net weight  | kg (lbs) |          | 170 (375)                     | 200 (441)   | 230 (507)   | 280 (617)   | 200 (441)                   | 280 (617)   |
| Shipping weight   | kg (lbs) |          | 190 (419)                     | 225 (496)   | 260 (573)   | 315 (694)   | 225 (496)                   | 315 (694)   |
| Max load capacity of one-piece work surface   | kg (lbs) |          | 50 (110)                      | 50 (110)    | 50 (110)    | 50 (110)    | 50 (110)                    | 50 (110)    |
| Max load capacity of single segmented work surfaces<br>(per 30cm segment)   | kg (lbs) |          | 25 (55)                       | 25 (55)     | 25 (55)     | 25 (55)     | 25 (55)                     | 25 (55)     |
| <b>Electrical Data</b>  |          |          |                               |             |             |             |                             |             |
| Voltage   | V        |          | 230                           |             |             |             | 120                         |             |
| Hertz   | Hz       |          | 50/60                         |             |             |             | 60                          |             |
| Protection class  |          |          | 1/IP20                        |             |             |             |                             |             |
| <b>Filter Specification</b>   |          |          |                               |             |             |             |                             |             |
| Specification   |          |          | H14 HEPA EN1822, 99.995% MPPS |             |             |             |                             |             |
| <b>Certifications</b>   |          |          |                               |             |             |             |                             |             |
| Standard  |          |          | EN 12469                      |             |             |             | NSF/ANSI 49                 |             |
| Agency certified  |          |          | TÜV *                         | TÜV         | TÜV *       | TÜV         | NSF                         |             |
| Certification LNE (France only)   | Option   |          |                               | LNE         |             |             | LNE                         |             |
| <b>Ergonomics and Utilities</b>   |          |          |                               |             |             |             |                             |             |
| Lighting intensity  | lux (fc) |          | 1000 (93)                     | 929 (86)    | 1500 (139)  | 1418 (132)  | 850 (79)                    | 1300 (121)  |
| Interior outlets  | pcs.     |          | 2 rear wall                   |             | 4 rear wall |             | 2 rear wall                 | 4 rear wall |
| Service valve ports   | pcs.     |          | 2 on each side + Smart Port   |             |             |             | 2 on each side + Smart Port |             |
| Slope of front window   | degree ° |          | 10                            |             |             |             | 10                          |             |
| <b>Herasafe 2030i Specifications</b>  |          |          |                               |             |             |             |                             |             |
| Noise level (measured according to EN 12469)<br>200mm opening   | dB (A)   |          | 57                            |             | 57          |             |                             |             |
| Noise level (EN units measured according to EN 12469,<br>NSF units measured according to NSF/ANSI 49)<br>250mm opening            |          |          | 57                            |             | 58          |             | 62                          | 62          |
| Power consumption (in normal operation, fans on,<br>lights on) 200mm opening  | W/h      |          | 133                           | 153         | 217         | 251         |                             |             |
| Power consumption (in normal operation,<br>Fans on, lights on) 250mm opening  |          |          | 155                           | 175         | 240         | 300         | 180                         | 316         |
| Power consumption (in stand-by mode)  | W/h      |          | 55                            | 55          | 80          | 80          | 55                          | 80          |
| Normal heat output in operation mode<br>(non-vented) 200mm opening  | W/h      |          | 133                           | 153         | 217         | 251         |                             |             |
| Normal heat output in operation mode<br>(non-vented) 250mm opening  |          |          | 155                           | 175         | 240         | 300         | 180                         | 316         |
| Stand-by heat output (lights off, non-vented)   | W/h      |          | 55                            | 55          | 80          | 80          | 55                          | 80          |

\* TÜV GS available for models with 200mm access opening

| Herasafe 2030i Ordering Information           |             |  | EN 12469  |           |           |           | NSF/ANSI 49 |            |
|---|-------------|--|-----------|-----------|-----------|-----------|-------------|------------|
| Size/Internal width                           | m (ft)      |  | 0.9 (3)   | 1.2 (4)   | 1.5 (5)   | 1.8 (6)   | 1.2 (4)     | 1.8 (6)    |
| Exhaust volume of cabinet 200mm opening       | m3/h (cf/m) |  | 296 (174) | 395 (232) | 493 (290) | 592 (348) |             |            |
| Exhaust volume of cabinet 250mm opening       | m3/h (cf/m) |  | 370 (218) | 494 (291) | 617 (363) | 741 (436) | 582 (342)   | 872 (513)  |
| Exhaust volume, thimble ducting 200mm opening | m3/h (cf/m) |  | 385 (226) | 513 (302) | 641 (377) | 770 (453) |             |            |
| Exhaust volume, thimble ducting 250mm opening | m3/h (cf/m) |  | 481 (283) | 642 (378) | 802 (472) | 963 (567) | 756 (445)   | 1134 (667) |

#### Receptacle Options

Country configurations are available to ensure correct line cord and internal receptacles are shipped with the cabinet. Standard version for EN12469 cabinets is Schuko, and for NSF cabinets is type B

|                               |        |  |          |  |  |  |                 |
|-------------------------------|--------|--|----------|--|--|--|-----------------|
| Country version Switzerland   | Option |  | 51900300 |  |  |  |                 |
| Country version UK            | Option |  | 51900303 |  |  |  |                 |
| Country version Italy         | Option |  | 51900306 |  |  |  |                 |
| Country version FR/B/CZ/SK/PL | Option |  | 51900921 |  |  |  |                 |
| Country version Australia     | Option |  | 51900449 |  |  |  |                 |
| Country version Denmark       | Option |  | 51900481 |  |  |  |                 |
| Country version China         | Option |  | 51900900 |  |  |  |                 |
| Country version Israel        | Option |  | 51901255 |  |  |  |                 |
| Country version India         | Option |  | 51901261 |  |  |  |                 |
| Country version Argentina     |        |  |          |  |  |  | 51900308 (230V) |
| Country version Chile         |        |  |          |  |  |  | 51901259 (230V) |
| Country version Brazil        |        |  |          |  |  |  | 51901260 (230V) |

#### Stand Options

We recommend ordering a base stand to ensure a robust and stable installation of the BSC. All of our base stands have been tested with our units for stability and functionality. Choose between a fixed height, manually adjustable working height, or electrically adjustable working height. Please consider the total height of the room when choosing the right working height. Footrest and castor options are available for added convenience.

|  |      |          |          |          |          |          |          |
|--|------|----------|----------|----------|----------|----------|----------|
| Fixed height stand charcoal grey<br>Height of 680 - 880 mm                           | Acc. | 50155665 | 50155666 | 50155667 | 50155668 | 50155666 | 50155668 |
| Adjustable height stand charcoal grey<br>Height of 680 - 880 mm                      | Acc. | 50155688 | 50155689 | 50155690 | 50155691 | 50155689 | 50155691 |
| Electric Adjustable Stand grey-white (RAL 9002)                                      | Acc. | 50117449 | 50109314 | 50117450 | 50109325 | 50109314 | 50109325 |
| Footrest bar   | Acc. | 50051982 | 50051983 | 50051984 | 50051985 | 50051983 | 50051985 |
| Castor-Set for Stand grey-white (4 pcs.) (RAL9002)<br>Provides work height of 850 mm | Acc. |          | 50109977 |          |          |          | 50109977 |

#### Work Surface Options

You have the option to fully customize your new BSC to accommodate any work task. The standard EN cabinet includes our stainless steel segmented worktray for easy removal, autoclaving, and cleaning under the work surface. The standard NSF cabinet includes our stainless steel single piece indented worktray. To customize your work area, order the part number to remove the standard worktray and configure the cabinet as desired. All internal surfaces and worktrays are 316 SS.

|  |              |          |          |          |          |          |          |
|--|--------------|----------|----------|----------|----------|----------|----------|
| Without table tops (negative price)                              | Option       | 51900651 | 51900652 | 51900653 | 51900654 | 51900652 | 51900654 |
| Flat worksurface module (150 mm)                                 | KAM15 Acc.   |          | 50073666 |          |          | 50073666 |          |
| Flat worksurface module (300 mm)                                 | KAM30 Acc.   |          | 50072685 |          |          | 50072685 |          |
| Flat worksurface module (600 mm)                                 | KAM60 Acc.   |          | 50073667 |          |          | 50073667 |          |
| Flat worksurface module (900 mm)                                 | KAM90 Acc.   |          | 50076398 |          |          | 50076398 |          |
| Lowered worksurface module (600 mm)                              | KAML60 Acc.  |          | 50073668 |          |          | 50073668 |          |
| Lowered worksurface module (900 mm)                              | KAML90 Acc.  |          | 50073669 |          |          | 50073669 |          |
| Lowered worksurface module (1200 mm)                             | KAML120 Acc. | No       | 50073670 |          |          | 50073670 |          |
| Lowered worksurface module (1500 mm)                             | KAML150 Acc. | No       | 50073671 |          |          | No       | 50073671 |
| Lowered worksurface module (1800 mm)                             | KAML180 Acc. | No       | 50078750 |          | No       | 50078750 |          |
| Vibration dampened surface module<br>with granite plate (300 mm) | KAMW30 Acc.  |          | 50073757 |          |          | 50073757 |          |
| Vibration dampened surface module<br>with granite plate (600 mm) | KAMW60 Acc.  |          | 50076162 |          |          | 50076162 |          |
| 3 liter mini-sink table top module                               | KAMS30 Acc.  |          | 50078875 |          |          | 50078875 |          |
| Replacement armrests KAS 30<br>(Set of 2 armrests 316 SS)        | KAS 30       |          | 50073663 |          |          | 50073663 |          |
| Padded armrest 316 SS<br>(Set of 2 including gel inserts)        | Acc.         |          | 50156438 |          |          | 50156438 |          |
| Replacement gel pads for padded armrest (2 pcs.)                 | Acc.         |          | 50156437 |          |          | 50156437 |          |

| Herasafe 2030i Ordering Information   |        |          | EN 12469                              |          |          |          | NSF/ANSI 49                           |          |
|---|--------|----------|---------------------------------------|----------|----------|----------|---------------------------------------|----------|
| Size/Internal width   | m (ft) |          | 0.9 (3)                               | 1.2 (4)  | 1.5 (5)  | 1.8 (6)  | 1.2 (4)                               | 1.8 (6)  |
| <b>Additional Filter Options</b>  |        |          |                                       |          |          |          |                                       |          |
| Additional exhaust HEPA filters (AEF) can be fitted to the top of the BSC to provide 99.9999975% cumulative exhaust filtration. Please note that AEF filters are not a requirement for recirculating BSC's, even in the UK, where BS5726 has been superceded by EN12469, since 2000. Additional exhaust carbon filters can be fitted on top of the BSC, but should only be used for traces of solvents being used as an adjunct to microbiological study. |        |          |                                       |          |          |          |                                       |          |
| Additional exhaust HEPA Filter, retro-fit   | KAEF   | Acc.     | 50078857                              | 50073772 |          | 50078856 |                                       |          |
| Additional exhaust Carbon Filter, retro-fit   | KAEFC  | Acc.     | 50078783                              | 50073773 |          | 50078784 |                                       |          |
| Grid filter set, G1 acc. To EN 779 (5 pcs)  |        |          | 50076846                              |          | 50076921 |          | 50076846                              | 50076921 |
| <b>UV Options</b>   |        |          |                                       |          |          |          |                                       |          |
| The UV control function is designed to preserve bulb life whilst monitoring constantly. The light is programmable via the user interface. This feature measures and tracks UV bulb use ensuring that the bulbs will not require excessive replacement due to premature and unnecessary degradation. Additionally by monitoring UV bulb age, the system auto-adjusts run time to ensure consistent sterilisation.  |        |          |                                       |          |          |          |                                       |          |
| Cross beam UV-C light option in side walls includes protective cover  |        |          |                                       |          |          |          |                                       |          |
| Included in models: 51032328, 51032329, 51032330, 51032331  |        | Standard | Standard (Cross Beam - in side walls) |          |          |          |                                       |          |
| Available for models: 51032324, 51032325, 51032326, 51032327  |        | Option   | 51900642 (Cross Beam - in side walls) |          |          |          |                                       |          |
| Included in models: 51032328, 51032329, 51032330, 51032331  |        | Standard |                                       |          |          |          | Standard (cross beam - in side walls) |          |
| Available for models: 51032324, 51032325, 51032326, 51032327  |        | Option   |                                       |          |          |          | 51900642 (cross beam - in side walls) |          |
| <b>Media Valves</b>   |        |          |                                       |          |          |          |                                       |          |
| Our standard BSCs feature a dual-side wall negative pressure design which allows the user to install service valves through the side access ports. We also offer factory-installed piping options on our premium units to connect utilities to the top of the BSC to conserve right to left exterior footprint.   |        |          |                                       |          |          |          |                                       |          |
| <b>Rear Wall Media Connections</b>  |        |          |                                       |          |          |          |                                       |          |
| 3 utility connectors prep (rear wall)<br>Select rear wall service taps additionally   |        | Option   |                                       |          |          |          |                                       |          |
| Available for models: 51032324, 51032325, 51032326, 51032327  |        | Option   | 51901251 (G3/8")                      |          |          |          |                                       |          |
| Included in models: 51032328, 51032329, 51032330, 51032331  |        | Standard |                                       |          |          |          | Included                              |          |
| Available for models: 51032324, 51032325, 51032326, 51032327  |        | Option   |                                       |          |          |          | 51901252 (NPT3/8")                    |          |
| <b>Rear Wall Service Taps/Media Valves (supplies are connected above BSC)</b>   |        |          |                                       |          |          |          |                                       |          |
| Combustible gas tap for rear wall   |        | Acc.     |                                       |          | 50068235 | 50068235 |                                       |          |
| Non-combustible gas tap for rear wall   |        | Acc.     |                                       |          | 50068234 | 50068234 |                                       |          |
| Water tap for rear wall   |        | Acc.     |                                       |          | 50068233 | 50068233 |                                       |          |
| Vacuum tap for rear wall  |        | Acc.     |                                       |          | 50068232 | 50068232 |                                       |          |
| <b>Side Wall Service Taps/Media Valves/Lead Throughs</b>  |        |          |                                       |          |          |          |                                       |          |
| Combustible gas tap for side wall   |        | Acc.     |                                       |          | 50046015 | 50046015 |                                       |          |
| Non-combustible gas tap for side wall   |        | Acc.     |                                       |          | 50045959 | 50045959 |                                       |          |
| Water tap for side wall   |        | Acc.     |                                       |          | 50044679 | 50044679 |                                       |          |
| Vacuum tap for side wall  |        | Acc.     |                                       |          | 50044678 | 50044678 |                                       |          |
| Hose connection for side wall, hose dia. 10-13mm  |        | Acc.     |                                       |          | 50059017 | 50059017 |                                       |          |
| Hose connection for side wall, hose dia. 3mm  |        | Acc.     |                                       |          | 50076408 | 50076408 |                                       |          |
| Hose connection for side wall, hose dia. 6/3mm  |        | Acc.     |                                       |          | 50066668 | 50066668 |                                       |          |
| <b>Solenoid Valve (should be quoted with combustible gas G3/8 connection)</b>   |        |          |                                       |          |          |          |                                       |          |
| Solenoid valve for combustible gas tap  |        | Acc.     |                                       |          | 50074905 | 50074905 |                                       |          |
| <b>Exhaust Duct Connection Kits</b>   |        |          |                                       |          |          |          |                                       |          |
| HEPA filters do a great job of trapping particles and biohazards, but they do not trap gases and toxic chemicals that may be used as an adjunct to microbiological studies. In these applications, it is important to select an exhaust duct transition to route exhaust air out of the laboratory to the external environment. We offer two different connections that can be installed in three different ways, depending upon facility requirements.   |        |          |                                       |          |          |          |                                       |          |
| Thimble Duct Exhaust Transition (Connection diameter 200mm/8")  |        | KTD      | Acc.                                  | 50078796 | 50073770 | 50078795 | 50073770                              | 50078795 |
| Direct Duct Exhaust Transition with anti blowback valve (Connection diameter 200mm/8")  |        | KDD      | Acc.                                  | 50078791 | 50073771 | 50078790 |                                       |          |
| Additional sensor and alarm to monitor flow in external duct. (Req. with thimble for NSF units)   |        |          |                                       | N/A      | N/A      | N/A      | 191085                                |          |



| Herasafe 2030i Ordering Information  |           |      | EN 12469 |          |         |          | NSF/ANSI 49 |         |
|--|-----------|------|----------|----------|---------|----------|-------------|---------|
| Size/Internal width  | m (ft)    |      | 0.9 (3)  | 1.2 (4)  | 1.5 (5) | 1.8 (6)  | 1.2 (4)     | 1.8 (6) |
| User Convenience Options   |           |      |          |          |         |          |             |         |
| We offer unique customer installed accessories to help organize your workspace and free up valuable space inside and around the BSC. |           |      |          |          |         |          |             |         |
| Hanging bar with set of 6 hooks  | IVB 6 6   | Acc. | 50073943 | No       |         |          | No          |         |
| Hanging bar with set of 12 hooks   | IVB 10 12 | Acc. | No       | 50073944 |         |          | 50073944    |         |
| Utility basket to be used with hanging bar   | KB        | Acc. | 50078900 |          |         | 50078900 |             |         |
| 2 utility baskets incl. additional hanging bar   | KBB       | Acc. | 50078901 |          |         | 50078901 |             |         |

| Herasafe 2030i Ordering Information  |  |  | Part number    |  | m (ft) |  |
|--|--|--|----------------|--|--------|--|
| Additional topboxes for fumigation   |  |  | Part number    |  | m (ft) |  |
| KSF frame fixed to top of each recirc cabinet to accept the relevant topbox for fumigation. Brackets are fitted to side of frame to allow the use of speed clamps. Height = 50mm   |  |  | KSF9           |  | Acc.   |  |
|  |  |  | KSF12/15       |  | Acc.   |  |
| KSNB topbox for use with Hydrogen Peroxide fumigation. Speed clamps allow for fixture to a KSF frame. The 35mm hole will accept a Steris OR Bioquell outlet camlock fitting. Height = 110mm with gasket.   |  |  | KSF18          |  | Acc.   |  |
|  |  |  | KSNB9          |  | Acc.   |  |
| KS-EVBOX for use with formonia system and formaldehyde. Speed clamps allow it to be clamped to KSF frame, while the 2 threaded holes accept the ABEKP3 charcol filters. This box can only be used if the formaldehyde is being neutralized with ammonia after fumigation. Height = 110mm with gasket.  |  |  | KSNB-12/15     |  | Acc.   |  |
|  |  |  | KSNB-18        |  | Acc.   |  |
| KSDB for use with formaldehyde fumigation only. Speed clamps allow for fixture to a KSP frame. This top box is supplied with a 3m flexi trunking which allows you to vent formaldehyde out into another exhausted application to purge fumes from cabinet. Side exit only. Height = 190mm with gasket. |  |  | KS-EVBOX-9     |  | Acc.   |  |
|  |  |  | KS-EVBOX-12/15 |  | Acc.   |  |
|  |  |  | KS-EVBOX-18    |  | Acc.   |  |
| KSP(FD) For use with formaldehyde fumigation. Fixes to a KSF frame. Supplied with a 3m flexi trunking which allows venting of fumes to an exhausted application. Top exit only. Height = 50mm (Frame).   |  |  | KS9DB160       |  | Acc.   |  |
|  |  |  | KS12DB160      |  | Acc.   |  |
|  |  |  | KS18DB160      |  | Acc.   |  |
| ABEKP3 charcoal filters screw into KS-EVBOX. Pack comprises of 3 filters.  |  |  | KSP(FD)-9      |  | Acc.   |  |
|  |  |  | KSP(FD)-12/15  |  | Acc.   |  |
|  |  |  | KSP(FD)-18     |  | Acc.   |  |
|  |  |  | ABEKP3         |  | Acc.   |  |

Not all BSCs are manufactured to the same safety standards. Our design philosophy is grounded in compliance with the recognized worldwide safety standards, because your work cannot compromise on safety. Our biological safety cabinets are subjected to all the following certifications:



EN 12469

Our Class II biological safety cabinets are independently tested and certified to EN 12469 standard for Class II biological safety cabinets by TÜV Nord, the leading testing agency in Europe.



LNE

Our Class II biological safety cabinets are independently tested and certified to EN 12469 by LNE, the leading testing agency in France.



NSF

Our Class II biological safety cabinets are independently tested and certified to NSF/ANSI 49, by NSF International, a leading testing agency in North America and recognized globally.

000 «Диаэм»

Москва  
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

**С.-Петербург**  
+7 (812) 372-6040  
spb@dia-m.ru

**Новосибирск**  
+7 (383) 328-0048  
nsk@dia-m.ru

**Воронеж**  
+7 (473) 232-4412  
vrn@dia-m.ru

**Йошкар-Ола**  
+7 (927) 880-3676  
nba@dia-m.ru

**Красноярск**  
+7 (923) 303-0152  
krsk@dia-m.ru

**Казань**  
+7 (843) 210-2080  
kazan@dia-m.ru

**Ростов-на-Дону**  
+7 (863) 303-5500  
rnd@dia-m.ru

**Екатеринбург**  
+7 (912) 658-7606  
ekb@dia-m.ru

**Кемерово**  
+7 (923) 158-6753  
kemerovo@dia-m.ru

**Армения**  
+7 (094) 01-0173  
armenia@dia-m.ru

