

Centrifuge compatibility for TaqMan Array Cards

| Current-generation centrifuges* | Cat. No. | | Cat. No. |
|--|------------|--|------------|
| Megafuge ST4 Plus-MD, 208–240 V, 50/60 Hz | 75009503 | Sorvall Legend XTR, 120 V (refrigerated) | 75004521 |
| Megafuge ST4 Plus, 208–240 V, 50/60 Hz | 75009903 | Sorvall Legend XTR, 230 V (refrigerated) | 75004520 |
| Megafuge ST4 Plus-MD, 120V | 75009504 | Sorvall Legend XTR, 230 V (US and Canada, refrigerated) | 75004523 |
| Megafuge ST4R Plus-MD, 220–240 V, 50 Hz/230 V, 60 Hz | 75009518 | Sorvall ST 40, 120 V | 75004510 |
| Megafuge ST4R Plus, 220 V 60 Hz | 75009818 | Sorvall ST 40, 230 V | 75004509 |
| Megafuge ST4R Plus-MD, 120 V | 75009519 | Sorvall ST 40R, 120 V (refrigerated) | 75004525 |
| Megafuge ST4R Plus, 220–240 V, 50 Hz/230 V, 60 Hz | 75009918 | Sorvall ST 40R, 230 V (refrigerated) | 75004524 |
| Megafuge ST4R Plus-MD, 220 V, 60 Hz | 75009618 | Sorvall ST4 Plus-MD, 208–240 V, 50/60 Hz | 75009509 |
| Multifuge X4 Pro, 208–240 V, 50/60 Hz | 75009900 | Sorvall ST4R Plus, 220 V, 60 Hz | 75009824 |
| Multifuge X4 Pro-MD, 120 V | 75009501 | Sorvall ST4R Plus-MD, 120 V | 75009525 |
| Multifuge X4 Pro-MD, 208–240 V, 50/60 Hz | 75009500 | Sorvall ST4R Plus-MD, 220 V, 60 Hz | 75009624 |
| Multifuge X4R Pro, 220–240 V, 50 Hz/230 V, 60 Hz | 75009915 | Sorvall X4 Pro, 208–240 V, 50/60 Hz | 75009905 |
| Multifuge X4R Pro-MD, 220–240 V, 50 Hz/230 V, 60 Hz | 75009515 | Sorvall X4 Pro-MD, 120 V | 75009506 |
| Multifuge X4R Pro, 220 V, 60 Hz | 75009815 | Sorvall X4 Pro-MD, 208–240 V, 50/60 Hz | 75009505 |
| Multifuge X4R Pro-MD, 120 V | 75009516 | Sorvall X4R Pro, 220–240 V, 50 Hz/230 V, 60 Hz | 75009920 |
| Multifuge X4R Pro-MD, 220 V, 60 Hz | 75009615 | Sorvall X4R Pro, 220 V, 60 Hz | 75009820 |
| Discontinued, previous-generation centrifuges* | Cat. No. | | Cat. No. |
| Megafuge 40, 120 V | 75004504 | SL 40R, 230 V (refrigerated) | 75004527 |
| Megafuge 40, 230 V | 75004503 | Sorvall Legend XT, 120 V | 75004506 |
| Megafuge 40R, 120 V (refrigerated) | 75004519 | Sorvall Legend XT, 230 V | 75004505 |
| Megafuge 40R, 230 V (refrigerated) | 75004518 | Sorvall Legend XT, 230 V (US and Canada) | 75004508 |
| Multifuge X3, 120 V | 75004501 | Sorvall ST4 Plus, 208–240 V, 50/60 Hz | 75009909 |
| Multifuge X3, 230 V | 75004500 | Sorvall ST4 Plus-MD, 120 V | 75009510 |
| Multifuge X3R, 120 V (refrigerated) | 75004516 | Sorvall ST4R Plus-MD, 220–240 V, 50 Hz/230 V, 60 Hz | 75009524 |
| Multifuge X3R, 230 V (refrigerated) | 75004515 | Sorvall ST4R Plus, 220–240 V, 50Hz/230 V, 60 Hz | 75009924 |
| SL 40, 120 V | 75004513 | Sorvall X4R Pro-MD, 120 V | 75009521 |
| SL 40, 230 V | 75004512 | Sorvall X4R Pro-MD, 220 V, 60 Hz | 75009620 |
| SL 40R, 120 V (refrigerated) | 75004528 | Sorvall X4R Pro-MD, 220–240 V, 50Hz/230 V, 60 Hz | 75009520 |
| Discontinued, previous-generation centrifuges** | Cat. No. | | Cat. No. |
| Multifuge 3S Plus, 120 V | 75004365 | Sorvall Legend T Plus, 120 V | 75004367 |
| Multifuge 3S Plus, 230 V | 75004361 | Sorvall Legend T Plus, 230 V | 75004363 |
| Multifuge 3SR Plus, 120 V (refrigerated) | 75004375 | Sorvall Legend RT Plus, 120 V | 75004377 |
| Multifuge 3SR Plus, 230 V (refrigerated) | 75004371 | Sorvall Legend RT Plus, 230 V (refrigerated) | 75004373 |
| Discontinued, previous-generation centrifuges† | Cat. No. | | Cat. No. |
| Fisher Scientific accuSpin Models 3/3R, 120 V | 13-100-517 | Fisher Scientific accuSpin Models 3/3R, 120 V (refrigerated) | 13-100-516 |

* Compatible with the Thermo Scientific™ TX-750 4 × 750 mL Swinging Bucket Rotor (75003180) and Applied Biosystems™ Bucket for TaqMan® Array Cards (4442571)

** Compatible with the Swing-Out Rotor TTH-750 4 × 750 mL (75006445) or Swing-Out Rotor 4 × 750 mL (75006434) and Applied Biosystems™ Bucket for TaqMan® Array Cards (4337762)

† Compatible with the centrifuge rotor (13-100-518) and Applied Biosystems™ Bucket for TaqMan® Array Cards (4337762)

000 «Диаэм»

Москва
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

