

Model: SUPREME LINE PLUS

Order number: 79-200960-021

Correct use: Compressed gas cylinders indoor storage



Classical door technology (door handle). For the storage of pressurised gas cylinders. Type testing by MPA/TÜV SÜD according to DIN EN 14470-2, DIN EN 16121 and DIN EN 16122, EK5/AK4 09-10:2009 and AfPS GS 2014:01 PAK Mark of approval: GS/CE mark, high quality mark of approval by TÜV SÜD Construction and colour:

- Outer carcass sheet steel, RAL 7035
- Doors sheet steel, RAL 3003
- Inner carcass of decor panels, RAL 7035 Fittings and features:
- Cylinder lock for safe access
- Roll-in flap with pitch
- Cylinder holder incl 3 retaining belts
- Fitting holder
- Exhaust air connection socket on the cabinet roof Further technical information: www.dueperthal.com

| | |
|-------------------|---------------------------|
| Category: | Type G90 - storage of gas |
| Model size | L |
| Outer width [mm]: | 895 |
| Outer depth [mm]: | 610 |

| | |
|---|----------------------|
| Outer height [mm]: | 2.070 |
| Net weight [kg]: | 480,00 |
| Number of storage levels | 1 |
| Number of storage areas per storage level | 1 |
| Door operation | classical technology |
| Door technology | wing door technology |

