

invitrogen

Streamline analysis, expand possibilities

E-Gel Power Snap Electrophoresis System

DNA separation

Simplify DNA electrophoresis with the only integrated gel running and imaging platform

The Invitrogen™ E-Gel™ Power Snap Electrophoresis System combines rapid, real-time nucleic acid analysis with high-resolution image capture for unmatched convenience. The integrated design helps reduce workflow time and accelerate discovery, offering:

- **Faster analysis**—go from sample loading to image capture in as little as 15 minutes
- **Simple operation**—intuitive user interface with large touch screen and integrated operating system
- **Safer handling**—minimize handling of hazardous chemicals when used with Invitrogen™ E-Gel™ precast gel cassettes



ДИА•М
современная лаборатория

www.dia-m.ru
заказ on-line

ThermoFisher
SCIENTIFIC

Benchtop convenience

The Invitrogen™ E-Gel™ Power Snap Electrophoresis Device is a small benchtop unit with a blue-light transilluminator for safer operation and improved downstream cloning efficiency. It also comes with an amber filter for real-time sample tracking on Invitrogen™ E-Gel™ agarose gels prestained with Invitrogen™ SYBR™ Safe or SYBR™ Gold II DNA stains. The device arrives with preprogrammed protocols for each type of E-Gel agarose gel.



For added convenience, stock your Supply Center with E-Gel cassettes and DNA ladders for direct and immediate access.

Request Invitrogen™ E-Gel™ products for your Supply Center at thermofisher.com/scproductrequests

Instant image capture

Capture high-resolution E-Gel agarose gel images quickly with the easy-to-use Invitrogen™ E-Gel™ Power Snap Camera. The camera docks directly to the electrophoresis unit, requiring no external power supply or connection to a desktop computer, enabling the freedom to capture images right at the bench.



Ordering information

| Product | Cat. No. |
|---|----------|
| E-Gel Power Snap Electrophoresis Device | G8100 |
| E-Gel Power Snap Electrophoresis Camera | G8200 |
| E-Gel Power Snap Electrophoresis System | G8300 |
| E-Gel Power Snap Electrophoresis System Starter Kit, EX 1% | G8341ST |
| E-Gel Power Snap Electrophoresis System Starter Kit, EX 2% | G8342ST |
| E-Gel Power Snap Electrophoresis Device Starter Kit, CloneWell II, 0.8% | G8168ST |
| E-Gel Power Snap Electrophoresis Device Starter Kit, SizeSelect II, 2% | G8162ST |

000 «Диаэм»

Москва
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

Новосибирск
пр. Академика
Лаврентьева, д. 6/1
тел.
(383) 328-0048
nsk@dia-m.ru

Казань
ул. Парижской
Коммуны, д. 6
тел.
(843) 210-2080
kazan@dia-m.ru

С.-Петербург
ул. Профессора
Попова, д. 23
тел.
(812) 372-6040
spb@dia-m.ru

**Ростов-
на-Дону**
пер. Семашко, д. 114
тел.
(863) 303-5500
rnd@dia-m.ru

Пермь
Представитель
тел.
(342) 202-2239
perm@dia-m.ru

Воронеж
Представитель
тел.
(473) 232-4412
voronezh@dia-m.ru

Армения
Представитель
тел.
(094) 01-0173
armenia@dia-m.ru

Узбекистан
Представитель
тел.
(90) 354-8569
uz@dia-m.ru