

(patented)


TRANSFER CAPS

**SIMPLE TRANSFER BETWEEN
FLASK, BAG & BIOREACTOR**

ДИАМ
современная лаборатория



www.dia-m.ru
заказ on-line


THOMSON™
Solutions
At Work
INSTRUMENT COMPANY



INVERSION FILL

BIDIRECTIONAL FILL







INVERSION FILL TRANSFER CAP

SAVE TIME & MAINTAIN
HIGH VIABILITY



**TRANSFER CAPS HAVE STANDARD
MALE QUICK CONNECT
FOR EASY CONNECTION TO BAGS**




**TRANSFER CAPS ARE AVAILABLE
WITH WELDABLE C-FLEX TUBING
(#16 OR #24) FOR STERILE
CONNECTION TO BAGS
OR BIOREACTORS**



**QUICK
CONNECT**









**TUBE
FUSING**




BIDIRECTIONAL PUMP TRANSFER CAP: TRANSFER CAPS WITH DOWN STEM FOR PRECISE CONTROL

The addition of the down stem tubing to the Optimum Growth Transfer Cap system allows for the Optimum Growth Flasks (patented) for use with a peristaltic pump. This design enables very precise control with the volume transferred to one or many vessels for inoculation from a single seed culture. The Pump Transfer Cap is bidirectional, so liquid can be pumped both in and out of the Optimum Growth flasks.



BIDIRECTIONAL CAPS



931614

Bidirectional Transfer Cap for Optimum Growth™ 5L Flask for Peristaltic Pump, Welds to 1/4" C-Flex 16 Tubing 2' -- Sterile Patented

Custom Tubing Available*


Please contact us at folks@htslabs.com

*Minimum order of 100 pieces required

Bidirectional Caps available in case quantities of 8


INVERSION CAPS

Quick Connect



931594


Inversion Transfer Cap for Optimum Growth™ 5L Flask with 1/4" Barb Connection -- Sterile



931596


Inversion Transfer Cap for Optimum Growth™ 5L Flask with 7/16" Male Connection -- Sterile

Tube Fusing



931595


Inversion Transfer Cap for Optimum Growth™ 5L Flask 1/4" Male Connection w/2' Tubing (Weld to C-Flex 16) -- Sterile



931598

Inversion Transfer Cap for Optimum Growth™ 5L Flask 7/16" Male Connection w/2' Tubing (Weld to C-Flex 24) -- Sterile

THOMSON
MADE IN USA



931607

Ring Only for Inverting Optimum Growth™ 5L Flask (must have 22" stand)



All Inversion Caps Require Ring & Stand

931606

Stand with Ring for Inverting Optimum Growth™ 5L Flask

000 «Диаэм»

Москва
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

