

ДИАМ  
современная лаборатория

www.dia-m.ru  
заказ on-line

eppendorf



# Get +Together

A perfect match: Centrifuge 5920 R + Rotor S-4xUniversal-Large + Eppendorf Consumables



**A perfect match**


Do you spin both conical tubes and plates? Eppendorf may have the perfect solution for you. Eppendorf offers a comprehensive range of high quality consumables and now introduces the new refrigerated Centrifuge 5920 R, featuring extraordinary high capacity and performance. Its main swing-bucket rotor allows for spinning both plates and tubes in the same bucket.

**Your benefits:**

- > Cost saving: No need for separate plate buckets
- > Space saving: No need to store different buckets
- > Time saving: No need to exchange buckets
- > Optimal complementarity and performance with Eppendorf Consumables

# The Power of +Together.

Find the perfect match for your application.



## Centrifuge 5920 R + Rotor S-4xUniversal-Large

- > Capacity up to 40 x 50 mL / 96 x 15 mL / 8 x DWP
- > Aerosol-tight QuickLock caps optionally available
- > Powerful state-of-the-art refrigeration system with advanced temperature management to keep your samples safe

## Eppendorf Tubes® + Eppendorf Plates®

- > High g-safe® centrifugation stability enables reliable, safe, and fast centrifugation steps
- > Different purity grades for individual application demands. Monitored and certified quality – lot by lot
- > Highest sample integrity due to avoidance of any slip agents, plasticizers, and biocides

# Get +Together online!



## Ordering information

Description	Max. capacity of conical tubes per rotor	International Order no.	North America Order no.
<b>Centrifuge 5920 R</b> , refrigerated	Depends on rotor	5948 000.018	5948000131
<b>Centrifuge 5920 R</b> , refrigerated, with Rotor S-4xUniversal-Large, incl. Universal buckets and adapters for 15 mL /50 mL conical tubes	96 x 15 mL/40 x 50 mL	5948 000.913	5948000930

Description	OptiTrack® frame color*	International Order no.	North America Order no.
<b>Eppendorf Conical Tubes</b> , 15 mL Sterile, pyrogen-, DNase-, RNase-, and DNA-free, 500 tubes (10 bags x 50 tubes)	—	0030 122.151	0030122151
<b>Eppendorf Conical Tubes</b> , 50 mL Sterile, pyrogen-, DNase-, RNase-, and DNA-free, 500 tubes (10 bags x 50 tubes)	—	0030 122.178	0030122178
<b>Microplate 96/F</b> , flat bottom, wells clear PCR clean, 80 plates (5 bags x 16 plates)	<input type="checkbox"/> white	0030 601.106	951040005
<b>Microplate 384/F</b> , flat bottom, wells clear PCR clean, 80 plates (5 bags x 16 plates)	<input type="checkbox"/> white	0030 621.107	951040341
<b>Microplate 96/F</b> , flat bottom, wells white PCR clean, 80 plates (5 bags x 16 plates)	<input checked="" type="checkbox"/> gray	0030 601.475	951040137
<b>Microplate 384/V</b> , V-shaped bottom, wells white PCR clean, 80 plates (5 bags x 16 plates)	<input checked="" type="checkbox"/> gray	0030 621.670	951040503
<b>Microplate 96/F</b> , flat bottom, wells black PCR clean, 80 plates (5 bags x 16 plates)	<input type="checkbox"/> white	0030 601.700	951040196
<b>Microplate 384/V</b> , V-shaped bottom, wells black PCR clean, 80 plates (5 bags x 16 plates)	<input type="checkbox"/> white	0030 621.905	951040481
<b>Deepwell Plate 96/500 µL</b> , wells clear, 500 µL PCR clean, 40 plates (5 bags x 8 plates)	<input type="checkbox"/> white	0030 501.101	951031801
<b>Deepwell Plate 96/1,000 µL</b> , wells clear, 1,000 µL PCR clean, 20 plates (5 bags x 4 plates)	<input type="checkbox"/> white	0030 501.209	951032603
<b>Deepwell Plate 96/2,000 µL</b> , wells clear, 2,000 µL PCR clean, 20 plates (5 bags x 4 plates)	<input type="checkbox"/> white	0030 501.306	951033405
<b>Deepwell Plate 384/200 µL</b> , wells clear, 200 µL PCR clean, 40 plates (5 bags x 8 plates)	<input type="checkbox"/> white	0030 521.102	951031003

\*More frame colors, purity grades, and bottom shapes are available at [www.eppendorf.com](http://www.eppendorf.com)

## 000 «Диаэм»

Москва  
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ [sales@dia-m.ru](mailto:sales@dia-m.ru)

[www.dia-m.ru](http://www.dia-m.ru)

**С.-Петербург**  
+7 (812) 372-6040  
[spb@dia-m.ru](mailto:spb@dia-m.ru)

**Новосибирск**  
+7 (383) 328-0048  
[nsk@dia-m.ru](mailto:nsk@dia-m.ru)

**Воронеж**  
+7 (473) 232-4412  
[vrn@dia-m.ru](mailto:vrn@dia-m.ru)

**Йошкар-Ола**  
+7 (927) 880-3676  
[nba@dia-m.ru](mailto:nba@dia-m.ru)

**Красноярск**  
+7 (923) 303-0152  
[krsk@dia-m.ru](mailto:krsk@dia-m.ru)

**Казань**  
+7 (843) 210-2080  
[kazan@dia-m.ru](mailto:kazan@dia-m.ru)

**Ростов-на-Дону**  
+7 (863) 303-5500  
[rnd@dia-m.ru](mailto:rnd@dia-m.ru)

**Екатеринбург**  
+7 (912) 658-7606  
[ekb@dia-m.ru](mailto:ekb@dia-m.ru)

**Кемерово**  
+7 (923) 158-6753  
[kemerovo@dia-m.ru](mailto:kemerovo@dia-m.ru)

**Армения**  
+7 (094) 01-0173  
[armenia@dia-m.ru](mailto:armenia@dia-m.ru)

