



Streamline your Cryogenic Storage

To maximize space and protect your precious samples

ΔΙΑ•M

Современная лаборатория

www.dia-m.ru

84922 00-0104

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand



Thermo Scientific™ Cryoboxes protect and organize your important samples and are available in materials ranging from economical cardboard to break-resistant, durable polycarbonate to meet any need and budget. The Dense Storage format is available to maximize storage, holding 169 1.0mL Thermo Scientific™ Nunc™ Cryobank™ tubes. All Thermo Scientific Cryoboxes are designed to provide dependable sample protection and meet the rigorous requirements of cryogenic storage.



Plastic Boxes

Made from break-resistant, durable polycarbonate, Thermo Scientific plastic cryoboxes protect samples during short- or long-term storage.

- Keyed lid prevents misalignment between lid and box
- Printed grid for easy sample management
- Autoclavable for decontamination
- Withstands temperature range of -196°C to +121°C
- Vented base ensures drainage and air circulation
- Available with dividers to hold 25 to 169 tubes
- Adaptable for cryogenic tubes ranging from 0.5mL to 5.0mL

Cardboard Boxes

Thermo Scientific white cardboard cryoboxes provide an economical choice for short-term sample storage.

- Available with dividers to hold 25 to 169 tubes
- Adaptable for cryogenic tubes ranging from 0.5mL to 4.5mL
- Allows handwriting with standard pens

| Cat. No. | Color box, lid | Tube Capacity, array | Tube Size, mL | Dimensions L x W x H mm (in) | Quantity per pack | Quantity per case |
|------------------------|-----------------|----------------------|---------------|-------------------------------------|-------------------|-------------------|
| White Plastic | | | | | | |
| 5025-0505 | White, Clear | 25, 5 x 5 | 1.0 – 2.0 | 76 x 76 x 52 (3.00 x 3.00 x 2.05) | 8 | 48 |
| 5026-0909 | White, Clear | 81, 9 x 9 | 1.2 – 2.0 | 132 x 132 x 52 (5.25 x 5.25 x 2.05) | 4 | 24 |
| 5026-1010* | White, Clear | 100, 10 x 10 | 1.2 – 2.0 | 132 x 132 x 52 (5.25 x 5.25 x 2.05) | 10 | 10 |
| 5027-0909 | White, Clear | 81, 9 x 9 | 5.0 | 132 x 132 x 95 (5.24 x 5.24 x 3.78) | 4 | 24 |
| 374187 | Gray, Gray | 100, 10 x 10 | 1.0 – 1.8 | 132 x 132 x 52 (5.25 x 5.25 x 2.05) | 4 | 24 |
| 341483 | Gray, Gray | 100, 10 x 10 | 3.6 | 132 x 132 x 78 (5.25 x 5.25 x 3.07) | 4 | 24 |
| Colored Plastic | | | | | | |
| 867013-0240 | Red, Clear | 81, 9 x 9 | 1.2 – 2.0 | 132 x 132 x 52 (5.25 x 5.25 x 2.05) | 4 | 24 |
| 867013-0241 | Yellow, Clear | 81, 9 x 9 | 1.2 – 2.0 | 132 x 132 x 52 (5.25 x 5.25 x 2.05) | 4 | 24 |
| 867013-0242 | Green, Clear | 81, 9 x 9 | 1.2 – 2.0 | 132 x 132 x 52 (5.25 x 5.25 x 2.05) | 4 | 24 |
| 867013-0243 | Blue, Clear | 81, 9 x 9 | 1.2 – 2.0 | 132 x 132 x 52 (5.25 x 5.25 x 2.05) | 4 | 24 |
| 867013-0244 | Gray, Clear | 81, 9 x 9 | 1.2 – 2.0 | 132 x 132 x 52 (5.25 x 5.25 x 2.05) | 4 | 24 |
| 867013-0245 | Assorted, Clear | 81, 9 x 9 | 1.2 – 2.0 | 132 x 132 x 52 (5.25 x 5.25 x 2.05) | 4 | 24 |

*Holds System 100 tubes

| Cat. No. | Color | Tube Capacity, array | Tube Size, mL | Dimensions L x W x H mm (in) | Quantity per pack | Quantity per case |
|--------------------------------|------------------|----------------------|---------------|-------------------------------------|-------------------|-------------------|
| Cardboard | | | | | | |
| 378247 | White | 25, 5 x 5 | 1.0 – 1.8 | 75 x 75 x 52 (2.88 x 2.88 x 2.00) | 24 | 24 |
| 378220 | White | 81, 9 x 9 | 1.0 – 1.8 | 132 x 132 x 52 (5.25 x 5.25 x 2.00) | 6 | 6 |
| 369639 | White | 81, 9 x 9 | 1.0 – 1.8 | 132 x 132 x 52 (5.25 x 5.25 x 2.00) | 48 | 48 |
| 369647 | White | 81, 9 x 9 | 3.6 | 132 x 132 x 75 (5.25 x 5.25 x 2.88) | 30 | 30 |
| 340053 | White | 100, 10 x 10 | 1.0 – 1.8 | 132 x 132 x 52 (5.25 x 5.25 x 2.00) | 42 | 42 |
| 340061 | White | 100, 10 x 10 | 3.6 | 132 x 132 x 75 (5.25 x 5.25 x 2.88) | 30 | 30 |
| Dense Storage Plastic | | | | | | |
| 5050-0169 | White, Clear lid | 169, 13 x 13 | 1.0 | 132 x 132 x 52 (5.25 x 5.25 x 2.00) | 4 | 24 |
| Dense Storage Cardboard | | | | | | |
| 369640 | White | 169, 13 x 13 | 1.0 | 132 x 132 x 52 (5.25 x 5.25 x 2.00) | 4 | 24 |

Compatibility Chart for Thermo Scientific Cryotubes and Boxes

| Cryotubes | | Plastic Cryoboxes | | | | | | Cardboard Cryoboxes | | | | |
|-------------------------|------------------------|-------------------|---|---------------------|---------------------|-----------|--------|---------------------|--------|------------------|--------|--|
| Catalog Number | Size, mL | 5025-0505 | 5026-0909 867013-0240 867013-0241 867013-0242 867013-0243 867013-0244 867013-0245 | 5026-1010 374187 | 341483 5027-0909 | 5050-0169 | 378247 | 378220 369639 | 369647 | 340053 340061 | 369640 | |
| Malgene Cryogenic Tubes | 5000-0012 5001-0012 | 1.2 | ● | ● | | | | ● | ● | | | |
| | 5000-0020 5001-0020 | 2.0 | ● | ● | | | | ● | ● | | | |
| | 5000-0050 5001-0050 | 5.0 | | | | ● | | | | | | |
| | 5000-1012 5001-1012 | 1.0 | ● | ● | ● | | | ● | ● | ● | | |
| | 5000-1020 5001-1020 | 1.5 | ● | ● | ● | | | ● | ● | ● | | |
| | 5011-0012 5011-0020 | 1.2 2.0 | ● ● | ● ● | | | | ● ● | ● ● | | | |
| | 5012-0012 5012-0020 | 1.2 2.0 | ● ● | ● ● | | | | ● ● | ● ● | | | |

| Cryotubes | | Plastic Cryoboxes | | | | | Cardboard Cryoboxes | | | | |
|----------------------|----------|-------------------|---|---------------------|-----------|--------|---------------------|--------|------------------|--------|--|
| Catalog Number | Size, mL | 5025-0505 | 5026-0909 867013-0240 867013-0241 867013-0242 867013-0243 867013-0244 867013-0245 | 5026-1010 374187 | 5050-0169 | 378247 | 378220 369639 | 369647 | 340053 340061 | 369640 | |
| Nunc Cryogenic Tubes | 374080 | 0.5 | ● | ● | | | ● | ● | | | |
| | 374115 | 1.0 | | | ● | | | | | ● | |
| | 375353 | 1.0 | ● | ● | | | ● | ● | | | |
| | 375299 | 1.0 | ● | ● | | | ● | ● | | | |
| | 347597 | 1.0 | ● | ● | | | ● | ● | | | |
| | 377224 | 1.0 | ● | ● | ● | | ● | ● | ● | | |
| | 366656 | 1.0 | ● | ● | ● | | ● | ● | ● | | |
| | 347651 | 1.0 | ● | ● | | | ● | ● | | | |
| | 347910 | 1.0 | ● | ● | | | ● | ● | | | |
| | 375418 | 1.8 | ● | ● | | | ● | ● | | | |
| | 340711 | 1.8 | ● | ● | | | ● | ● | | | |
| | 347627 | 1.8 | ● | ● | | | ● | ● | | | |
| | 377267 | 1.8 | ● | ● | ● | | ● | ● | ● | | |
| 368632 | 1.8 | ● | ● | | | ● | ● | | ● | | |
| 363401 | 1.8 | ● | ● | ● | | ● | ● | | ● | | |

Compatibility Chart for Thermo Scientific Cryotubes and Boxes

| Cryotubes | | Plastic Cryoboxes | | | | | | | Cardboard Cryoboxes | | | | |
|--------------------------------|----------|-------------------|---|--|---------------------|--------|-----------|-----------|---------------------|------------------|--------|------------------|--------|
| Catalog Number | Size, mL | 5025-0505 | 5026-0909 867013-0240 867013-0241 | 867013-0242 867013-0243 867013-0244 867013-0245 | 5026-1010 374187 | 341483 | 5027-0909 | 5050-0169 | 378247 | 378220 369639 | 369647 | 340053 340061 | 369640 |
| Nunc Cryogenic Tubes Continued | 351934 | 1.8 | ● | ● | ● | | | | ● | ● | | ● | |
| | 347783 | 1.8 | ● | ● | ● | | | | ● | ● | | ● | |
| | 379189 | 3.6 | | | | ● | ● | | | | ● | | |
| | 366524 | 3.6 | | | | ● | ● | | | | ● | | |
| | 367997 | 3.6 | | | | ● | ● | | | | ● | | |
| | 337516 | 4.5 | | | | | ● | | | | | | |
| | 347643 | 4.5 | | | | | ● | | | | | | |
| | 379146 | 4.5 | | | | | ● | | | | | | |
| | 363452 | 4.5 | | | | | ● | | | | | | |
| 347791 | 4.5 | | | | | ● | | | | | | | |
| Nunc Universal Tubes | 374502 | 1.8 | ● | ● | | | | | ● | ● | | | |
| | 374503 | 1.8 | ● | ● | | | | | ● | ● | | | |
| | 374512 | 2.0 | ● | ● | | | | | ● | ● | | | |
| | 374513 | 2.0 | ● | ● | | | | | ● | ● | | | |

Compatibility Chart for Thermo Scientific Cryoboxes and Freezer Racks

| Cryoboxes Racks | Vertical Racks | | | | | | Horizontal Racks | | Adjustable | |
|-----------------|----------------|-----------------------------|-----------------------------|---------------------|---------------------|-----------------------|-----------------------|----------------------|------------|----------------------|
| | Catalog Number | DS5035-0004 (5025-0505 x 4) | DS5035-0009 (5025-0505 x 9) | 5036-0004 (2mL x 4) | 5036-0009 (2mL x 9) | DS5037-0004 (5mL x 4) | DS5037-0007 (5mL x 7) | 5038-4322 (2mL x 12) | | 5038-4422 (2mL x 16) |
| 5025-0505 | ● | ● | | | | | | | | |
| 5026-0909 | | | ● | ● | | | | ● | ● | ● |
| 5026-1010 | | | ● | ● | | | | ● | ● | ● |
| 5027-0909 | | | | | ● | ● | | | | |
| 374187 | | | ● | ● | | | | ● | ● | ● |
| 341483 | | | ● | ● | | | | ● | ● | ● |
| 867013-0240 | | | ● | ● | | | | ● | ● | ● |
| 867013-0241 | | | ● | ● | | | | ● | ● | ● |
| 867013-0242 | | | ● | ● | | | | ● | ● | ● |
| 867013-0243 | | | ● | ● | | | | ● | ● | ● |
| 867013-0244 | | | ● | ● | | | | ● | ● | ● |
| 867013-0245 | | | ● | ● | | | | ● | ● | ● |
| 378247 | ● | ● | | | | | | | | |
| 378220 | | | ● | ● | | | | ● | ● | ● |
| 369639 | | | ● | ● | | | | ● | ● | ● |
| 369647 | | | | | ● | ● | | | | |
| 340053 | | | ● | ● | | | | ● | ● | ● |
| 340061 | | | ● | ● | | | | ● | ● | ● |
| 5050-0169 | | | ● | ● | | | | ● | ● | ● |
| 369640 | | | ● | ● | | | | ● | ● | ● |



Peace of Mind with Each Lot

Thermo Scientific cryoboxes are backed by the Thermo Fisher Scientific high quality production standards, customer service, and technical support that make Thermo Scientific the most trusted name in cryogenic storage.

Visit our cryogenic storage resource page at [thermoscientific.com/cryotubes](https://www.thermoscientific.com/cryotubes) to learn more about our sample storage solutions or contact your sales representative.



000 «Диаэм»

Москва

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург

+7 (812) 372-6040
spb@dia-m.ru

Новосибирск

+7(383) 328-0048
nsk@dia-m.ru

Воронеж

+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола

+7 (927) 880-3676
nba@dia-m.ru

Красноярск

+7(923) 303-0152
krsk@dia-m.ru

Казань

+7(843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону

+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург

+7 (912) 658-7606
ekb@dia-m.ru

Кемерово

+7 (923) 158-6753
kemerovo@dia-m.ru

Армения

+7 (094) 01-0173
armenia@dia-m.ru

